


# StentSpect<sup>CR</sup>




When it comes to Inspection, we have the Vision

## ➔ Key Features


**F**ully automatic non contact, accurate and repeatable crimped stent imaging and visual documentation system




**D**isplay stitched image of entire crimped stent (or balloon)



**G**raphical presentation of profile diameter measurement



**R**efine manufacturing process through graphical SPC analytical report



## ➔ Specifications

- ➔ Measures crimped stent up to 55 mm Length (other options available)
- ➔ Measures crimped stent up to 5 mm diameter (other options available)
- ➔ Scans average long length crimped stent in less than 30 seconds
- ➔ Accuracy: Better than 6 micron

Distributed by:

